



VOXELLO[®]
BETTER OUTCOMES.

Voxello Training

Supporting Patient-Provider Communication Improves Patient Outcomes

Patients who face communication barriers are at increased risk of preventable adverse events, for example:

- Falls
- Pressure Ulcers
- Ventilator Associated Pneumonias
- Adverse Drug or Blood Product Reactions

To make patients active participants in their care.

- Patients need to be able to summon their nurses
 - Access the nurse call system
- Patients need to be able to effectively communicate with their nurses
 - To communicate about symptoms
 - To participate in assessments
 - To indicate their preferences about cares
 - To ask questions about cares

Barriers to Patient-Provider Communication

- Inability to use conventional nurse call pendant
- Inability to use alternative nurse call system switches (e.g. soft touch)
- Inability to speak due to intubation or tracheostomy
- Inability to produce intelligible speech
- Inability to write due to weakness or paresis

Using Assistive Technology to Overcome Barriers.

- Take advantage of any intentional gesture a patient can make
- Enable patients to access the nurse call system with a small intentional gesture
- Enable patients to use intentional gestures to control a speech generating device

The Voxello System




The First Part Of Patient-Provider Communication: Accessing The Nurse Call System

- All conscious patients should have access to the nurse call system.
 - Standard nurse call pendant
 - Standard alternative switches
 - Pressure plate
 - Pressure bulb
- Patients who lack the motor skills to use the call pendant (pillow speaker) or alternative switches need some form of assistive technology to access the call system.

BETTER OUTCOMES.



Nurse Call Decision Tree



- Identify the voluntary gesture that the patient can produce.
 - What is a gesture?
 - Any motor response a patient can make on command
 - Any motor response a patient can make repeatedly
 - Types of gestures
 - tongue click
 - pushing tongue into cheek
 - lip squeeze
 - wink or eye blink
 - Head, jaw, shoulder or limb movement

Selecting a transducer to detect the intentional gesture


- Gestures with an audible component
 - Microphone
- Gestures with a minimal movement and very limited force
 - Proximity sensor
- Gestures with minimal movement and force
 - Minmo (gyro/accelerometer) sensor or Pressure sensor
- Winking or blinking gestures
 - Infrared wink sensor

The Second Part Of Patient-Provider Communication: Accessing A Speech-Generating Device

- All conscious patients should have a way to communicate with their caregivers and participate in their care.
- Speech-generating devices (SGD) provide an option for patients who cannot use normal oral and written forms of communication.
 - **Direct Selection** using a touch screen for patients with use of their hands
 - **Scan Selection** using switches for patients with limited motor skills.

Communication Decision Tree

- Some patients, due to their medical status, may benefit from one sensor at one point in time and a different sensor at a later point in time.
- It is important to continually check on the patient to ensure a good fit between a patient's abilities and the selected sensor.



Who can benefit from using Voxello's technology

- Patients above the age of 3
- Patients who are conscious and not continuously sedated
- Patients who can produce an intentional gesture
- Patients who can understand that producing an intentional gesture
 - can activate the nurse call system
 - control a speech generating device

Noddle® Setup & Instructions



noddle Powered On



noddle Charging Port



noddle Input Indicating
Sensor Not Plugged In



noddle Input Sensor
Indicating Gesture
Detected

Using the noddle® to access the nurse call



noddle Sensors & Mounting


- J Touch
- Bed Touch
- VL Mic
- J Mic
- *Minmo*
- *Wink*



J Mic Sensor & J Touch Sensor

- The J Mic Sensor & J Touch Sensor both have cervical collar and head strap mounting capabilities through a Velcro strap.

- When using the noddle-touch position the sensor near to the patient's cheek and ensure that it can be activated when the patient's tongue is pressed against his or her cheek.



Bed Touch Sensor

- The noddle Bed Touch can be clipped to bedding with clips or mounted on a pillow using a band.
- In both cases, care must be taken to ensure that the patient is able to access the sensor.



VL Mic Sensor

- The noddle VL Mic can only be attached to the vent line for patients who are endotracheally intubated.
- The VL mic should easily snap on the vent line. Ensure that the microphone is properly positioned so that the patient can easily activate the noddle with a tongue click.



Using the noddle To Access The noddle-chat





- The noddle connected to noddle-chat via Bluetooth®.
 - → *tablet settings* → *Bluetooth* → *scan* → *paired devices* → *choose noddle*
- Use two-switch row-column scanning to navigate.
- Use gestures to scan menu
 - 1 to step through options
 - 2 to select to select an option

noddle-chat Setup

- New User: Select the appropriate user (child, teen, adult & gender) from the list when setting up a new user.

- Vocabulary Set Options:
Vocabulary sets include Start, Basic, Full + and Full - .



Kiosk Mode

'Kiosk Mode' locks the user in noddle-chat, so that the user can not navigate out of the program.

To set Kiosk mode:

- Choose 'Settings'
- Choose 'System'
- Choose 'Kiosk Mode'
- Check 'Kiosk Mode'



Voxello Bedside Setup





www.voxello.com





VOXELLO[®]
BETTER OUTCOMES.